
QIAGEN Patents

Matching "QIAGEN" by Owner.

IPMonitor
Patents

Contents

Alerts	3
"QIAGEN"	3
Terms and Conditions	6
General	6
Disclaimer of warranty and limitation of liability	6
Copyright	6
Arbitration	6

Alerts

"QIAGEN"

83 results matching "QIAGEN" by Owner.

Number	Title	Owner
2002314455	Preparation of polymer particles	Qiagen GmbH
2002314457	Preparation of polymer particles	Qiagen GmbH
2002315302	DOUBLE-STRAND NUCLEIC ACIDS WITH ONE OR TWO INDIVIDUAL STRAND 5' OVERHANGS, METHOD FOR PRODUCTION AND USE THEREOF	QIAGEN GMBH
2002245018	Method and device for collecting and stabilizing a biological sample	BECTON DICKINSON AND COMPANY ; QIAGEN GMBH
2002247654	Ammonium sulfate for neutralization of inhibitory effects	Qiagen GmbH
2002346960	METHODS AND COMPOSITIONS FOR PROMOTING GROWTH AND INNATE IMMUNITY IN YOUNG ANIM- ALS	QIAGEN GMBH ; UNIVERSITY OF SASKATCHEWAN
2002360930	METHOD FOR ISOLATING DNA FROM BIOLOGICAL SAMPLES	QIAGEN GMBH
2002211719	Compositions and methods for using a solid support to purify RNA	Qiagen North American Holdings, Inc.
2003208740	Method for coarse purification of cell digests from microorganisms	Qiagen GmbH
2003210362	METHOD FOR MODIFYING THE CONCENTRATION OF TRANSCRIPTS IN BIOLOGICAL SAMPLES CONTAINING RIBONUCLEIC ACID	QIAGEN GMBH
2003216572	Process for isolating nucleic acid with chaotrope agents and ammonium compounds	QIAGEN AS
2003219396	ISOLATING NUCLEIC ACID	QIAGEN AS
2003230070	ISOLATING NUCLEIC ACID	QIAGEN AS
2004202359	Mixing and pouring apparatus	Qiagen North American Holdings, Inc.
2001263800	Porous ferroor ferrimagnetic glass particles for isolating molecules	Qiagen GmbH
2001271722	Signal amplification with lollipop probes	Yale University ; Qiagen, GmbH
2003304421	Sample presentation device with differing wettability	Qiagen Sciences, Inc.
2004221678	Method and kit for a specific detection of M.tuberculosis	Qiagen Diagnostics GmbH
2004233035	Compositions and methods for using a solid support to purify RNA	Qiagen North American Holdings, Inc.
2004243576	Method for obtaining plasmid-DNA by means of an aqueous biphasic system	Qiagen GmbH
2004259851	Method for the reverse transcription and/or amplification of nucleic acids	Qiagen GmbH
2001268725	Multiply-primed amplification of nucleic acid sequences	Qiagen GmbH
2005209389	Method for chromatographic separation of a nucleic acid mixture	Qiagen GmbH

Number	Title	Owner
2005246390	Sample presentation device	Qiagen Sciences, Inc.
2001269031	Novel compositions for isolating and/or stabilising nucleic acids in biological material	Qiagen GmbH
2005268888	Monolithic shaped bodies for purifying and separating biopolymers	Qiagen GmbH
2007200429	Compositions and methods for using a lysing matrix for isolating DNA	Qiagen North American Holdings, Inc.
2007201583	Novel compositions for isolating and/or stabilising nucleic acids in biological material	Qiagen GmbH
2005305012	Compositions and methods for purifying nucleic acids from stabilization reagents	Qiagen North American Holdings, Inc.
2007203372	Universal collection medium	QIAGEN Gaithersburg, Inc.
2006212473	Sample lysis and coating of reaction surface	Qiagen GmbH
2006212392	Method for isolating nucleic acids, the nucleic acids being immobilised on a matrix at an increased temperature	Qiagen GmbH
2007216758	Mixing and pouring apparatus	Qiagen North American Holdings, Inc.
2006248997	Extraction of proteins from formalin-fixed tissue	Qiagen GmbH
2006298738	Method for processing a fluid and fluid processing device	Qiagen GmbH
2006298739	Apparatus for processing biological material	Qiagen GmbH
2006290435	Modulation of immunostimulatory properties of short interfering ribonucleic acid (siRNA) by nucleotide modification	Coley Pharmaceutical GmbH ; Qiagen GmbH
2006311288	Detection of nucleic acids by target-specific hybrid capture method	QIAGEN Gaithersburg, Inc.
2006323937	Magnetic polymer particles	Qiagen GmbH
2006323954	Method for enriching short-chain nucleic acids	Qiagen GmbH
2006334404	A method for treating a biological sample	Qiagen GmbH
2007217549	Method for production of a transformed cell	Qiagen GmbH
2007278425	Device for processing samples	Qiagen GmbH
2007286300	Method for synthesis of a cDNA in a sample in an enzymatic reaction	Qiagen GmbH
2007341280	Method and materials for triggered release of a biological sample	Qiagen GmbH
2008220769	Fixation of a biological material	QIAGEN GmbH
2008235605	Method for purifying biomolecules	Qiagen GmbH
2008288398	Method for suspending or resuspending particles in a solution and apparatus adapted thereto	Qiagen GmbH
2008290533	Apparatus for processing biological material	Qiagen GmbH
2008290599	Apparatus and process for removing substances from pre-filled containers	Qiagen GmbH
2008290600	Apparatus for delivering pipettable substances	Qiagen GmbH
2008281720	Method and device for fixing/stabilising a sample	Qiagen GmbH
2008286254	Method for synthesizing a cDNA in a sample in an enzymatic reaction	Qiagen GmbH
2009238247	Compositions, methods, and kits using synthetic probes for determining the presence of a target nucleic acid	QIAGEN Gaithersburg, Inc. ; Nazarenko, I. ; O'Neil, D.
2009237945	Fluorescence standard, and the use thereof	QIAGEN Lake Constance GmbH

Number	Title	Owner
2009243574	Simultaneous detection of multiple nucleic acid sequences in a reaction	QIAGEN Hamburg GmbH ; QIAGEN GmbH
2009254307	Method for isolating nucleic acids	QIAGEN GmbH
2009253197	Lysis, binding and/or wash reagent for isolating and/or purifying nucleic acids.	QIAGEN GmbH
2009257644	Magnetic microplate assembly	QIAGEN Gaithersburg, Inc.
2009300084	Method of loading a crystallization device	QIAGEN GmbH ; Forschungszentrum Jülich GmbH(FJZ)
2009294755	Method and apparatus for the simultaneous, automated decomposition of a plurality of biological samples	QIAGEN GmbH
2009294655	Method for normalizing the contents of biomolecules in a sample	QIAGEN GmbH
2009290982	Reaction vessel for crystallizing a sample from a solution	Qiagen GmbH
2009306339	Quantification of RNA using internal normalization	QIAGEN GmbH
2011201979	Processes for isolating, amplifying and characterizing DNA	Qiagen North American Holdings, Inc
2009320213	Fast results hybrid capture assay and system	QIAGEN Gaithersburg Inc.
2009332887	Nucleic acid purification method	QIAGEN GmbH
2009332900	Nucleic acid purification method	QIAGEN GmbH
2009353641	Closure and method of using same	QIAGEN
2010213126	RNAi based selection system	QIAGEN GmbH
2010208311	Sequence-specific large volume sample preparation method and assay	QIAGEN Gaithersburg, Inc.
2010208386	Thermophilic helicase dependent amplification technology with endpoint homogenous fluorescent detection	QIAGEN Gaithersburg, Inc.
2010244453	Detection of multiple nucleic acid sequences in a reaction cartridge	Qiagen GmbH
2010242867	A non-target amplification method for detection of RNA splice-forms in a sample	QIAGEN Gaithersburg, Inc.
2010246587	Nucleic acid purification method	QIAGEN GmbH
2009347207	Ensuring sample adequacy using turbidity light scattering techniques	QIAGEN
2010298620	Compositions, methods, and kits for isolating and analyzing nucleic acids using an anion exchange material	QIAGEN Gaithersburg Inc. ; QIAGEN GmbH
2010291990	Compositions and methods for recovery of nucleic acids or proteins from tissue samples fixed in cytology media	QIAGEN Gaithersburg, Inc.
2011203450	Methods, compositions, and kits for recovery of nucleic acids or proteins from fixed tissue samples	QIAGEN Gaithersburg, Inc.
2011204313	Materials and methods for isothermal nucleic acid amplification	QIAGEN Gaithersburg, Inc.
2012205204	Method for enriching short-chain nucleic acids	Qiagen GmbH
2011210734	Methods and compositions for sequence-specific purification and multiplex analysis of nucleic acids	QIAGEN Gaithersburg, Inc.
2011210748	Method of determining and confirming the presence of HPV in a sample	QIAGEN Gaithersburg, Inc.

Terms and Conditions

General

The information provided in this report is not in the nature of legal or other professional advice.

The information provided in this report is derived from third parties and may contain errors. You must make your own enquiries and seek independent advice from the relevant industry professionals regarding any information in this report.

Patents, patent applications, trademarks and trademark applications shown as lapsed or ceased may be restored at a later date.

Disclaimer of warranty and limitation of liability

Neither IPMonitor, its affiliates nor any of their respective employees, agents, third party content providers or licensors warrant that the information contained in this report is error free; nor do they make any warranty as to the accuracy, reliability or content of any information in this report.

This report is provided on an "as is" basis; to the maximum extent permitted by law, IPMonitor disclaims all representations and warranties, express or implied, with respect to any information contained in this report. In addition, IPMonitor does not represent or warrant that the information contained in this report is accurate, complete or current.

This disclaimer of liability applies to any damages or injury caused by any failure of performance, error, omission, interruption, deletion, defect, delay in operation or transmission, computer virus, communication line failure, theft or destruction or unauthorized access to, alteration of, or use of record, whether for breach of contract, tortious behavior, negligence, or under any other cause of action.

Copyright

The content of this report may be protected by copyright laws and treaties around the world and all rights relating to the same are reserved. Any paper or digital copies of this report are not to be modified or any illustrations or any graphics used separately from any accompanying text.

Arbitration

Any dispute or difference whatsoever arising out of or in connection with this report shall be submitted to arbitration in accordance with, and subject to, The Institute of Arbitrators & Mediators Australia Rules for the Conduct of Commercial Arbitrations.

IPMonitor reserves the right to update the Terms and Conditions at any time without notice. Please see the full Terms and Conditions at www.IPMonitor.com.au/trademarks/info/terms